

IL-12 p35 + p70 Antibody : B-T21 Azide free

NEUTRALISING ANTI HUMAN INTERLEUKIN-12 p35 + p70 MONOCLONAL ANTIBODY

CATALOGUE N° : 855.120.005

CLONE : B-T21

ISOTYPE : Mouse IgG1

IMMUNISATION : Recombinant human IL-12

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises both natural and recombinant human IL-12 p35 + p70

QUANTITY : 0.5 mg / 0.5 ml

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Store at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Blocks IL-12 induced proliferation on PHA and IL-2 activated lymphocytes
Suitable for intra-cytoplasmic flow cytometry

REFERENCES : Nagai, T. et al., J Immunol.,2003; 171(10): 5233-43.
Oppmann, B. et al.,Immunity,2000; 13(5): 715-25.

Version 2 - 01.08

For research use only