

**IL-12 p40 + 70 Antibody : B-P40 Azide free**

NON NEUTRALISING ANTI HUMAN INTERLEUKIN-12 p40 + p70 MONOCLONAL ANTIBODY

**CATALOGUE N° :** 855.180.005

**CLONE :** B-P40

**ISOTYPE :** Mouse IgG1

**IMMUNISATION :** Recombinant human IL-12

**HYBRIDOMA :** Myeloma X63/AG.8653 x Balb/c spleen cells

**SPECIFICITY :** Recognises both natural and recombinant human IL-12 p40 + p70

**QUANTITY :** 0.5 mg / 0.5 ml

**PURIFICATION :** Ion exchange chromatography

**FORMAT :** Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

**STORAGE :** Store at +2-8°C for 12 months  
For longer storage, freeze aliquots

**APPLICATION :** No effect on IL-12 induced proliferation on PHA and IL-2 activated lymphocytes  
Suitable for intra-cytoplasmic flow cytometry

**REFERENCES :** Verreck, F.A.W. et al., PNAS 2004; 101 : 4560 - 4565  
**Version 3- 04.07** **For research use only**