

IL-13 Antibody : B-B13 Azide free

NEUTRALISING ANTI HUMAN INTERLEUKIN-13 MONOCLONAL ANTIBODY

CATALOGUE N° : 855.140.005

CLONE : B-B13

ISOTYPE : Mouse IgG1

IMMUNISATION : Recombinant human IL-13

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognises both natural and recombinant human IL-13

QUANTITY : 0.5 mg / 0.5 ml

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Store at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Blocks IL-13 induced proliferation on TF1 cell line
Suitable for intra-cytoplasmic flow cytometry

REFERENCES : Guillot C. et al., Blood, 2003; 101 : 3325 - 3333
van der Pouw Kraan, T. C. et al., J Immunol.,1996; 156(5):
1818-23.

Woerly, G. et al.,J Leukoc Biol.,2002; 72(4): 769-79.

Version 2 - 01.08

For research use only