

### IL-13 Antibody : B-B13 PE conjugated

ANTI HUMAN INTERLEUKIN-13 PE CONJUGATED MONOCLONAL ANTIBODY

**CATALOGUE N° :** 855.142.010

**CLONE :** B-B13

**ISOTYPE :** Mouse IgG1

**IMMUNISATION :** Recombinant human IL-13

**HYBRIDOMA :** Myeloma X63/AG.8653 x Balb/c spleen cells

**SPECIFICITY :** Recognises both natural and recombinant IL-13

**QUANTITY :** 100 µg/1ml


**PURIFICATION :** Ion exchange chromatography

**FORMAT :** Phosphate-buffered saline with 1% BSA and 0,1% sodium azide

**WORKING DILUTION :** Reconstitute with 1 ml deionized water  
It is recommended to establish the optimal dilution according to the system being used. Use 10µl of antibody to label  $5 \cdot 10^5$  cells

**STORAGE :** Store at +2-8°C for six months after reconstitution  
DO NOT FREEZE

**APPLICATION :** Suitable for intra-cytoplasmic flow cytometry

 *Sodium azide in a poisonous and hazardous substance, handle and dispose of properly*

**Version 4 - 10.07**

**For research use only**