

IL-17A Antibody : B-B51 Azide free

NEUTRALISING ANTI HUMAN IL-17A AZIDE FREE

CATALOGUE N° : 853.910.000

CLONE : B-B51

ISOTYPE : Mouse IgG2b

IMMUNISATION : Recombinant human IL-17A

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c node cells

SPECIFICITY : Recognizes both recombinant and natural human IL-17A

QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Neutralises the bioactivity of IL-17A on keratinocytes cell line HaCaT (measurement of β defensin-2 production)